

TIMS

High precision tungsten isotopes analyzed via thermal ionization mass spectrometry

Authors

Andrea Mundi-Petermeier¹,
Jenny Roberts², Hauke Vollstaedt²

¹University of Vienna, Austria

²Thermo Fisher Scientific, Bremen,
Germany

Keywords

Tungsten isotopes, TIMS

Goal

To demonstrate the high precision achievable for tungsten isotope measurements using the negative mode of the Thermo Scientific™ Triton™ XT thermal ionization mass spectrometer.

Introduction

The short-lived ^{182}Hf - ^{182}W isotopic system [$^{182}\text{Hf} \rightarrow ^{182}\text{W} + 2\beta^-$, half-life ($t_{1/2}$) = 8.9 million years]¹ is a powerful tool for studying planetary differentiation processes that occurred within the first 60 million years of solar system history. Given that tungsten (W) is moderately siderophile and hafnium (Hf) is strongly lithophile, the system is particularly sensitive to metal-silicate segregation. The most widely used application for this system has been dating core formation of planetary bodies.²⁻⁴ By measuring W isotopic compositions of iron meteorites, ages can be calculated, assuming chondritic Hf/W. The ^{182}Hf - ^{182}W isotopic system has also been used for terrestrial applications. $^{182}\text{W}/^{184}\text{W}$ is used for measuring ^{182}W ingrowth from lithophile Hf. For example, the W isotopic compositions of mantle-derived rocks have been used to investigate mantle evolution.⁵⁻⁸ Isotopic variations in $^{182}\text{W}/^{184}\text{W}$ can be smaller than 10 ppm,^{6,9} requiring high precision W isotope measurements.

Here we demonstrate the capabilities of the Thermo Scientific Triton XT thermal ionization mass spectrometer equipped with $10^{13} \Omega$ Amplifier Technology to measure W isotopes at high precision of better than 2 ppm (2 RSD) for 0.3–1.25 μg W loadings.

Experimental

Sample preparation

Single Re filaments (99.99% purity) were used in these measurements. Filaments were outgassed at 4.5 A for at least 48 hours prior to loading, ideally several weeks.

Filament loading procedures generally followed Archer et al.¹⁰ Specifically, W was loaded in HCl matrix on the filament followed by a brief heating of the filament to a dull glow. To enhance ionization, 1 µL of activator solution containing 20 µg of La and 5 µg of Gd in 1 M Teflon double-distilled HCl was then added in 2–3 aliquots to the standard.

Instrumentation

Analyses were performed at the University of Vienna on a Triton XT system in negative ion mode. A combination of 10¹¹ and 10¹³ Ω amplifiers were used during analysis (Table 1). Specifically, the low noise 10¹³ Ω amplifiers were used to measure the small ion beams of ¹⁸⁰W¹⁶O₃, ¹⁸¹Ta¹⁶O₃, ¹⁸⁶W¹⁶O₂¹⁸O, ¹⁸⁷Re¹⁶O₂¹⁸O, and ¹⁹⁰Os¹⁶O₃ during the multi-static method.

A total of 18 blocks were measured with 20 cycles per block for a total measurement time of ~11 hours. The six 10¹¹ Ω amplifiers were rotated after every block (a total of three full amplifier rotations across the run). Peak centering was performed

every 6th block, a lens focusing every 2nd block and a 6 minute baseline was performed every block. High-precision 160 min amplifier gains were determined before every other sample run; however, in the data reduction, an average of the last three gain measurements was used to correct each 10¹¹ Ω amplifier.

Oxide production was enhanced using a Varian™ leak valve to bleed oxygen (P_{source} = 1.2 × 10⁻⁷ mbar) into the source⁹.

Data reduction

Data reduction procedures follow a three-step iterative process¹⁰ including interference correction of Re and Ta isotopologues and oxygen corrections. Potential Os interferences were also monitored via the measurement of ¹⁹⁰Os¹⁶O₃ using 10¹³ Ω amplifiers on H4 in line 2 (Table 1). All ratios were corrected for instrumental mass fractionation using a ¹⁸⁶W/¹⁸³W ratio of 1.9859. A two-sigma outlier rejection was made based on cycles.

Results and discussion

The internal precision of the individual measurements is between 2.1 ppm and 2.6 ppm (2 RSE) and close to the limit of counting statistics (Table 2). The results of two independent analytical sessions have ¹⁸²W/¹⁸⁴W external precisions of 2.1 and 1.8 ppm (2 RSD; n = 7 for each session; Table 2; Figure 1).

Table 1. Cup configuration for W isotope measurements

10 ¹³ Ω	10 ¹¹ Ω	10 ¹¹ Ω	10 ¹¹ Ω	10 ¹¹ Ω	10 ¹¹ Ω	10 ¹¹ Ω	10 ¹³ Ω	10 ¹³ Ω	Integration time	Idle time
L4	L3	L2	L1	CC	H1	H2	H3	H4	[s]	[s]
¹⁸⁰ W ¹⁶ O ₃	¹⁸¹ Ta ¹⁶ O ₃	¹⁸² W ¹⁶ O ₃	¹⁸³ W ¹⁶ O ₃	¹⁸⁴ W ¹⁶ O ₃	¹⁸⁵ Re ¹⁶ O ₃	¹⁸⁶ W ¹⁶ O ₃	¹⁸⁶ W ¹⁶ O ₂ ¹⁸ O	¹⁸⁷ Re ¹⁶ O ₂ ¹⁸ O	33.554	10
¹⁸¹ Ta ¹⁶ O ₃	¹⁸² W ¹⁶ O ₃	¹⁸³ W ¹⁶ O ₃	¹⁸⁴ W ¹⁶ O ₃	¹⁸⁵ Re ¹⁶ O ₃	¹⁸⁶ W ¹⁶ O ₃	¹⁸⁷ Re ¹⁶ O ₃	¹⁸⁷ Re ¹⁶ O ₂ ¹⁸ O	¹⁹⁰ Os ¹⁶ O ₃	33.554	10

Table 2. ¹⁸²W/¹⁸⁴W ratios for 14 analyses from two analytical sessions of a 0.3-1.25 µg Alfa Aesar™ W standard

Wheel	Sample	Average ¹⁸⁴ W signal [V]	¹⁸² W/ ¹⁸⁴ W multi static	2σ	2 RSE (ppm)	Theoretical limit of precision 2 RSE (ppm)
1	1	1.4	0.864861	2.10E-06	2.4	2.0
1	2	1.2	0.864861	2.27E-06	2.6	2.2
1	3	1.3	0.864863	2.47E-06	2.9	2.1
1	4	1.3	0.864860	2.23E-06	2.6	2.1
1	5	1.1	0.864862	2.07E-06	2.4	2.3
1	6	1.2	0.864861	2.14E-06	2.5	2.2
1	7	1.1	0.864860	2.25E-06	2.6	2.2
		Average	0.864861	1.78E-06	2.1	
2	1	1.2	0.864861	2.19E-06	2.5	2.2
2	2	1.1	0.864860	2.60E-06	3.0	2.2
2	3	1.2	0.864860	2.37E-06	2.7	2.2
2	4	1.2	0.864860	2.11E-06	2.4	2.2
2	5	1.2	0.864858	2.11E-06	2.4	2.2
2	6	1.2	0.864860	2.27E-06	2.6	2.2
2	7	1.1	0.864860	2.23E-06	2.6	2.3
		Average	0.864860	1.59E-06	1.8	




Figure 1. $^{182}\text{W}/^{184}\text{W}$ ratios for 14 analyses from two analytical sessions of a 0.3–1.25 μg Alfa Aesar W standard. Mean $^{182}\text{W}/^{184}\text{W}$ and 2 SD for each session shown by solid and dashed lines. Error bars represent 2 SE internal precisions.

Based on the average signal intensity of each filament (given a measurement time of 402 min = 2 lines x 33.554 s * 18 blocks * 20 cycles), the theoretical limit on internal precision of the measurement can be determined:

$$2\text{RSE (ppm)} = \sqrt{\left(\frac{2 \times 10^6}{\sqrt{\text{Tot. counts}_{^{184}\text{W}}}}\right)^2 + \left(\frac{2 \times 10^6}{\sqrt{\text{Tot. counts}_{^{182}\text{W}}}}\right)^2}$$

Table 2 shows that the measurements made on the 14 W samples were within 40–50% of the theoretical limit of internal precision based on counting statistics. This demonstrates that the precision of the method is not limited by the instrument noise but rather by signal intensity and total measurement time.

Conclusion

The Triton XT TIMS system equipped with $10^{13} \Omega$ Amplifier Technology is capable of analyzing $^{182}\text{W}/^{184}\text{W}$ isotope ratios at precisions of ≤ 2 ppm for 0.3–1.25 μg W loadings. This allows small W isotopic variations to be identified, opening up new possibilities for exploring small-scale heterogeneities resulting from early planetary differentiation processes. These isotopic differences can help in constraining the timing of mantle homogenization in Early Earth or tracing core-mantle interactions, ultimately giving a better understanding of geodynamic processes and composition of mantle reservoirs.

Acknowledgments

A M-P acknowledges support from FWF grant V659-N29.

References

- Vockenhuber, C. et al. New half-life measurement of ^{182}Hf : improved chronometer for the early solar system. *Phys Rev Lett* 2004, 93, 172501.
- Kruijer, T. S. et al. Protracted core formation and rapid accretion of protoplanets. *Science* 2014, 344, 1150–1154.
- Harper, C. L. & Jacobsen, S. B. Evidence for ^{182}Hf in the early Solar System and constraints on the timescale of terrestrial accretion and core formation. *Geochimica et Cosmochimica Acta* 1996, 60, 1131–1153.
- Lee, D. C. & Halliday, A. N. Hafnium–tungsten chronometry and the timing of terrestrial core formation. *Nature* 1995, 378, 771–774.
- Mundl-Petermeier, A. et al. Anomalous ^{182}W in high $3\text{He}/4\text{He}$ ocean island basalts: Fingerprints of Earth's core? *Geochimica et Cosmochimica Acta* 2020, 271, 194–211.
- Mundl, A. et al. Tungsten-182 heterogeneity in modern ocean island basalts. *Science* 2017, 356, 66–69.
- Willbold, M.; Mojzsis, S. J.; Chen, H. W. & Elliott, T. Tungsten isotope composition of the Acasta Gneiss Complex. *Earth and Planetary Science Letters* 2015, 419, 168–177.
- Willbold, M.; Elliott, T. & Moorbath, S. The tungsten isotopic composition of the Earth's mantle before the terminal bombardment. *Nature* 2011, 477, 195–198.
- Mundl-Petermeier, A.; Viehmann, S.; Tusch, J.; Bau, M. & Münker, C. Earth's geodynamic evolution constrained by ^{182}W in Archean seawater. *Nature Communications*, in press, 2022, doi: 10.1038/s41467-022-30423-3.
- Archer, G. J.; Mundl, A.; Walker, R. J.; Worsham, E. A. & Bermingham, K. R. High-precision analysis of $^{182}\text{W}/^{184}\text{W}$ and $^{183}\text{W}/^{184}\text{W}$ by negative thermal ionization mass spectrometry: Per-integration oxide corrections using measured $^{180}/^{160}$. *Int J Mass Spectrom* 2017, 414, 80–86.

Learn more at thermofisher.com/tims